CONDUCTIMETRO PORTÁTIL PORTABLE CONDUCTIVITY METER

Modelo / Model 910/3

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This manual should be available for all users of these equipments. To get the best results and a higher duration of this equipment it is advisable to read carefully this manual and follow the processes of use.

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ANNEX I: CE CERTIFICATE

AUXILAB S.L.

CE DECLARATION OF CONFORMITY
NAHITA PORTABLE CONDUCTIVITY METER 910/3 of AUXILAB, S.L. for the
Directive of Machines (89/392/CEE modified)
and the regulations adopted for their transposition

NAME OF THE MANUFACTURER / IMPORTER:

AUXILAB, S.L.

ADDRESS:

Polígono Morea Norte, 8
31191 Beriáin (Navarra)

WE STATE THAT:

NAHITA PORTABLE CONDUCTIVITY METER 910/3
CODE 50910030

Is designed and manufactured according to:

- ♦ Directive 89/392/CEE, including the modifications and national regulations that transpose them.
- ♦ Directive 89/336/CEE modified over electromagnetic compatibility.

And that the following harmonized rules have been applied (or part of them):

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BERIAIN 19th June 2007

Signed by: ALFONSO AINCIBURU SANZ
 DIRECTOR/MANAGER

Polígono Morea Norte, 8 31191 Beriain (Navarra) - Spain. Tel. 948 310 513 Fax 948 312 071
Internet: www.auxilab.es · Email: correo@auxilab.es

Cleaning

- ◆ Never use scourers or substances that can grate for cleaning the case of the equipment as they damage the conductivity meter and produce an early ageing of the equipment.
- ◆ Use a fluff-free cloth dampened with soaped water that does not contain abrasives.

ATTENTION!

IF EQUIPMENTS ARE NOT PROPERLY CLEAN AND DISINFECTED THEY WOULD NOT BE ALLOWED TO REPAIR BY OUR TECHNICAL SERVICE.

INSTRUCTIONS ON ENVIRONMENT PROTECTION

At the end of its life cycle, please, do not dispose of this equipment by throwing it in the usual garbage; hand it over a collection point for the recycling of electrical and electronic appliances. It does not contain dangerous or toxic products for humans but a non adequate disposal would damage the environment.

The materials are recyclable as mentioned in its marking. By recycling material or by other forms of re-utilization of old appliances, you are making an important contribution to protect our environment.

Please inquire at the community administration for the authorized disposal location.

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
Security

- ◆ The conductivity meter must be used by previously qualified staff that knows how the equipment works thanks to the user manual.
- ◆ It is essential to have the equipment switched off before cleaning, checking components or replacing any piece.
- ◆ Never try to repair the conductivity meter by yourself, since you will lose the warranty and may provoke damages to the general operating system or the electrical installation, as well as injuries to the people that usually handle the equipment (burns, hurts...).
- ◆ Try not to spill any liquid on the control panel, though it is properly insulated.
- ◆ Made under the European regulations for electrical security, electromagnetic compatibility and security on machines.

Accessories

- ◆ Nahita conductivity meter 910/3 is supplied with a conductivity cell with a constant of 1 cm^{-1} . In case of cell breakage or if a conductivity cell with a different constant value is needed, the following models are available in our catalogue:

Code	Constant	Range
90908010	$0.1 \pm 0.02 \text{ cm}^{-1}$	0.1 μS -300 mS/cm
90908020	$1 \pm 0.2 \text{ cm}^{-1}$	1 μS -100 mS/cm
90908030	$10 \pm 2 \text{ cm}^{-1}$	300 μS -1000 mS/cm

5. MAINTENANCE AND CLEANING

To get the best results and a higher duration of the conductivity meter it is essential to follow the processes of use.

Note: All the processes of use mentioned below will not have any value unless you keep a continued and careful maintenance.

- ◆ Please follow the processes of use of this manual.
- ◆ This manual should be available for all users of this equipment.
- ◆ Always use original components and supplies. Other devices can be similar but they can damage the equipment.
- ◆ Never try to repair the conductivity meter by yourself, since you will lose the warranty and may provoke damages to the general operating system, as well as injuries to the people that usually handle the conductivity meter (burns, hurts...) or damages in nearby equipments.
- ◆ Maintenance of the conductivity cell. Rinse the conductivity cell in clean water after each measurement; for a thorough cleaning, clean the cell with a cloth or a non abrasive detergent. After cleaning the cell, calibrate the equipment.
- ◆ Do not use the light more than necessary to extend the life of batteries.
- ◆ In the event of breakdown please contact your distributor to overhaul through Nahita Technical Assistance Department.

Calibration of the conductivity cell constant

- ◆ Each conductivity cell presents a specific constant value; if user thinks that the value indicated by the factory is not correct, proceed as follows:
- ◆ Select the right standard solution depending on the constant value of the cell:

KCl standard solution to be used depending on the cell constant value

Cell constant	0.01	0.1	1	10
Approx. concentration of KCl solution (mol/L)	0.001	0.001	0.01 or 0.1	0.1 or 1

Chart for preparing the different KCl standard solutions

Concentration (mol/L)	KCL concentration (g/L) at 20°C
1	74.2650
0.1	7.4365
0.01	0.7440
0.001	Dissolve 100 mL of the 0.01 mol/L solution in 1 L

Note: The KCl solution must be heated at 110 °C for 4 h, and then let cool.

- ◆ Connect the conductivity cell to the equipment and the conductivity meter on.
- ◆ Clean the conductivity cell with distilled water.
- ◆ Submerge the conductivity cell in the standard solution. Control the temperature of the standard solution to be (25.0±0.1) °C.
- ◆ Press ▲(4) or ▼(5) to adjust the temperature value on display to 25 °C. Press Enter (3) to save the value.
- ◆ Press Mode (2) and with ▲ (4) and ▼(5) adjust the constant value until the conductivity reading matches the value indicated at the following chart according to the standard solution used and its temperature.

Temperature	Approximate concentration (mol/L)			
	1	0.1	0.01	0.001
	Conductivity S/cm			
15°C	0.09212	0.010455	0.0011414	0.0001185
18°C	0.09780	0.011163	0.0012200	0.0001267
20°C	0.10170	0.011644	0.0012737	0.0001322
25°C	0.11131	0.012852	0.0014083	0.0001465
35°C	0.13110	0.015351	0.0016876	0.0001765

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The following chart shows the recommended cell constant value according to the measuring range:

Conductivity range ($\mu\text{S/cm}$)	Resistivity range (Ω/cm)	Recommended conductivity cell constant (cm^{-1})	Modification of the results to obtain the real value of conductivity
0-0.2 0-2 0-20 0-200	∞ - 5000000 ∞ - 500000 ∞ - 50000 ∞ - 5000	0.01	Displayed result x 0.01
0-20 0-200	∞ - 50000 ∞ - 5000	1	Displayed result x 1
0-2000 0-10000	∞ - 500 ∞ - 100	1	Displayed result x 1
0-200 0-2000 0-20000 0-100000	∞ - 5000 ∞ - 500 ∞ - 50 ∞ - 10	10	Displayed result x 10

Working temperature adjustment

- ♦ Nahita conductivity meter 910/3 is provided with a temperature compensation mechanism by entering data using the keypad. For this purpose:
- ♦ In the measuring mode press ▲ (4) or ▼ (5) for a few seconds until the °C symbol starts blinking.
- ♦ Enter the exact value of the temperature of the solution to be analyzed by using ▲ (4) and ▼ (5).
- ♦ Press *Enter* (3) to save the temperature value.

Conductivity measurement

- ♦ Turn the conductivity meter on by pressing *On/Off* (1).
- ♦ Clean the conductivity cell with distilled water and dry it with a filter or cellulose paper.
- ♦ Submerge the cell and the temperature probe 4 cm into the sample to be analyzed and stir it with the cell. Wait until the reading is stable. Conductivity is displayed at the middle and below is shown the temperature value previously entered.

Note: If several samples are going to be analyzed successively it is recommended to clean the conductivity cell and the temperature probe with distilled water to avoid cross contamination.

PLEASE READ THOROUGHLY THE INSTRUCTIONS BEFORE OPERATING WITH THIS EQUIPMENT WITH THE AIM OF ACHIEVING THE HIGHEST FEATURES AND THE MAXIMUM DURABILITY OF THE INSTRUMENT.

♦ Nahita conductivity meter 910/3 operates with 2 alkaline batteries AAtype. Open the batteries compartment located at the back side of the equipment and put the batteries according to their polarity. Put again the lid of the compartment.

Note: Remove the batteries if you are not using the instrument for a long period of time.

♦ Connection of the conductivity cell. Connect the conductivity cell(8) to the connector (7) placed at the upper part of the equipment.

♦ Turn the conductivity meter on by pressing *On/Off* (1). After a few seconds the equipment will be ready to operate.

Note: In order to assure precise results, it is necessary to adjust and establish some parameters (see *Parameter adjustment*) before performing any reading.

Parameter adjustment

Conductivity cell constant adjustment

♦ The instrument is supplied with a conductivity cell with a constant value of 1 cm^{-1} . To adjust the exact value of the cell constant (supplied by the manufacturer), proceed as following:

♦ Press *Mode* (2).

♦ With \blacktriangle (4) and \blacktriangledown (5) adjust the value that is shown on display to the exact value of the constant (E.g.: if the exact value of the constant is 0.95, adjust the displayed value to .950).

♦ Press *Enter* (3) to save the value. The equipment will automatically return to the measuring mode.

IMPORTANT: The conductivity meter can also be used with conductivity cells with a constant value of 0.01 or 10 cm^{-1} (see *Accessories*). To adjust the exact constant value of this kind of conductivity cells, proceed as follows:

♦ Press *Mode* (2).

♦ With \blacktriangle (4) and \blacktriangledown (5) adjust the value shown on display to the exact value of the cell constant.

♦ Cells with a constant value of 0.01 cm^{-1} . If the exact value of the constant is 0.012 adjust the value displayed to 1.200. This way, the result of the measurements should be multiplied by 0.01 to obtain the real conductivity value.

♦ Cells with a constant value of 10 cm^{-1} . If the exact value of the constant is 11 adjust the value displayed to 1.100. This way, the result of the measurements should be multiplied by 10 to obtain the real conductivity value.

♦ Press *Enter* (3) to save the value. The equipment will automatically return to the measuring mode.

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3. TECHNICAL SPECIFICATIONS

Code	50910030
Conductivity:	
Range	0.00-19.99 µS/cm; 20.0-199.9 µS/cm; 200-1999 µS/cm; 2000-19000 µS/cm; 20.0-100.0 mS/cm
Resolution	0.01 µS/cm; 0.1 µS/cm; 1 µS/cm; 0.1 mS/cm
Accuracy (± 1 digit)	± 1.5%
Temperature range	15-35 °C
Temperature accuracy	± 1.0%
Temperature compensation	Entering data using the keypad
Temperature coefficient	2 %/°C (fixed)
Reference temperature	25 °C
Working conditions	5-35 °C, = 85% HR
Power	2 batteries 1.5 V, AA type
Dimensions (LxAxH)	170x80x35 mm
Weight	0.3 Kg

4. INSTALLATION / SETTING UP**Preliminary inspection**

- ♦ Unwrap the conductivity meter and make sure that the equipment does not present any damage because of the shipment. In case it presents any damage, tell it immediately to your transport agent or dealer so that they can make the claims in the correct time limit.
- ♦ Please keep the original wrapping; you will always need it for returns enclosed with all the accessories supplied.
- ♦ Please check that all the accessories are enclosed with the equipment:
 - Conductivity cell
 - 2 batteries 1.5 V, AA type
 - User's manual
 - Warranty certificate

We will only accept any equipment return within 15 days after delivery and provided it comes in its original wrapping.

Installation

- ♦ Before using the conductivity meter, it is convenient for you to familiarize with its components and basic essentials, as well as with its control functions.

INDEX OF CONTENTS


1. USES OF THE INSTRUMENT.....	12
2. DESCRIPTION	12
3. TECHNICAL SPECIFICATIONS.....	13
4. INSTALLATION / SETTING UP	13
5. MAINTENANCE AND CLEANING	17
ANNEX I: CE CERTIFICATE	19

1. USES OF THE INSTRUMENT

Nahita conductivity meter 910/3 is a light, compact and easy-to-use portable conductivity meter ideal for *in situ* conductivity determinations. It presents an easy-to-read liquid crystal display that can be illuminated to facilitate readings in dark places and a damp-proofed membrane keypad with 5 functional keys to access to all the instrument possibilities.

The equipment is supplied with a cell with two platinum measuring plates ($K=1 \text{ cm}^{-1}$), and allows measuring conductivity of solutions in a wide range from 0 to $1 \times 10^5 \text{ } \mu\text{S/cm}$. Moreover, the value of the cell constant can be modified to allow the use of other conductivity cells with a constant value of 10, 1 and 0.01 cm^{-1} , thus covering a large variety of applications. It also presents a temperature compensation mechanism by entering data using the keypad.

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Thank you for choosing this equipment. We sincerely wish that you enjoy your Nahita portable conductivity meter 910/3. We highly recommend looking after this equipment according to what is stated in this manual.

Nahita develops its products according to the CE marking regulations as well as emphasizing the ergonomics and security for its user. The correct using of the equipment and its good quality will permit you to enjoy this equipment for years.

The improper use of the equipment can cause accidents and electric discharges, circuit breakers, fires, damages, etc. Please read the point of *Maintenance*, where we expose the security notes.

TO GET THE BEST RESULTS AND A HIGHER DURATION OF THE EQUIPMENT IT IS ADVISABLE TO READ THOROUGHLY THIS MANUAL BEFORE OPERATING WITH THE EQUIPMENT.

Please bear in mind the following:

- ◆ This manual is inseparable from the Nahita portable conductivity meter 910/3, so it should be available for all the users of this equipment.
- ◆ You should carefully handle the equipment avoiding sudden movements, knocks, free fall of heavy / sharp objects on it. Avoid spilling liquids inside the equipment.
- ◆ Never dismantle the different pieces of the equipment to repair it yourself, since it could produce a defective use of the whole equipment and a loss of the product warranty, as well as injuries on people that handle the conductivity meter.
- ◆ To prevent fire or electric discharges avoid dry or dusty environments. In case it may happen unplug the equipment immediately.
- ◆ If you have any doubt about setting up, installation or functioning do not hesitate in contacting your wholesaler. You can also tell us any doubts or suggestions you have by contacting Nahita Technical Assistance Department by email to asistencia@auxilab.es or by telephone: +34 807 117 040 (0.30 Euros/min).
- ◆ This equipment is protected under the Warranties and consumer goods regulation (10/2003).
- ◆ Overhaul is not covered by the equipment warranty.
- ◆ Operations made by non-qualified staff will automatically produce a loss of the warranty.
- ◆ Accessories, including their loss are not covered by the product's warranty. The warranty neither covers piece's deterioration due to the course of time.
- ◆ Please make sure you keep the invoice, either for having the right to claim or asking for warranty coverage. In case you have to send the equipment to Nahita Technical Assistance Department you should enclose the original invoice or a copy as guarantee.
- ◆ Manufacturer reserves the right to modify or improve the manual or equipment.

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